

**MASSACHUSETTS INSTITUTE OF TECHNOLOGY**  
Department of Electrical Engineering and Computer Science

**6.161 Modern Optics Project Laboratory**

Laboratory Exercise No. 1

Fall 2011

**Safety, Polarization States of Optical Waves, Dielectric Reflection and Geometric Optics**

---

To get the most out of your in-lab experience, you must come to Lab prepared (makes life easier for you and the TA, and minimizes your time in the Lab). Thus, you should go through this Lab manual, complete the Pre-Lab Exercises, and answer all the Pre-Lab questions **BEFORE** entering the Laboratory. In your lab notebook record data, explain phenomena you observe, and answer the questions asked. Remember to answer all questions in your lab notebook in a neat and orderly fashion. No data are to be taken on these laboratory sheets. Tables provided herein are simply examples of how to record data into your laboratory notebooks. Expect the in-lab portion of this exercise to take about 3 hours. **Instead of a formal written report, your report for this Laboratory Exercise will be presented in oral fashion before the TA and the Writing Coordinator. Detailed instructions are supplied in Appendix I.** The grading rubric is shown in Appendix II.

---

## PRE-LAB EXERCISES

### PL1.1 - Safety

Read the entire Laboratory and laser Safety Packet that was handed out in Lab 0. Be prepared to answer questions about high-voltage and laser safety before arriving at the MOL.

### PL1.2 – – Get Prepared to Start the Laboratory Exercises

Read the **entire** laboratory handout, and be prepared to answer questions before, during and after the lab session. Determine all the equations and constants that may be needed in order to perform all the laboratory exercises. **Write** them all down in your laboratory notebook before entering the Lab. This will ensure that you take all necessary data while in the Lab in order to complete the lab write-up. This preparatory work will also count toward your Lab Exercise grade.

### PL1.3 – Telescope Design

For this exercise, you are encouraged to work in groups of 2-3 people. You are to design a telescope given the optical components listed in Table I. Your telescope must be a refractor (lenses-only). The entire telescope assembly must not be any longer than the 2 feet (the shorter, the better). Your telescope will need to contain **at least** 2 optical elements, but no more than 4. The telescope must be a terrestrial telescope (meaning that the image as seen by a viewer must be erect). Your total system magnification must be at least 4x and at most 40x at an imaging distance of either 15 meters or adjusted to infinity (a collimating telescope – this is the most common form of telescope, and is preferred for simplicity and robustness of design).

Be sure to use the ABCD matrix approach to simplify the design of your telescope. You should include all of your system specifications in the pre-lab write-up – including any computer code, hand-calculations or other materials used in the design. Be sure to specify the focal lengths, total system magnification as well as the locations of the lenses, the object and the image. Remember, there is no one ‘correct’ way to design your telescope. The Lab staff will be very sad if everyone designs the same telescope. Additional points will be awarded for ingenuity and thoughtfulness in design. Also, be sure not to use too many elements; each additional element in your system will lead to decreased contrast and decreased resolution. “Ingenuity and thoughtfulness” means designing smart, and minimizing complexity.

TABLE I. Optical elements available for building your telescope:

Size	Type	Focal Length	Quantity
2"	Bi-Convex	200mm	1
2"	Bi-Convex	100mm	2
2"	Bi-Convex	75mm	1
1"	Bi-Convex	200 mm	1
1"	Bi-Convex	150 mm	1
1"	Bi-Convex	100 mm	1
1"	Bi-Convex	75 mm	1
1"	Bi-Convex	50 mm	1
1"	Bi-Convex	35 mm	2
1"	Bi-Convex	25.4 mm	2
1"	Bi-Concave	-100mm	1
1"	Bi-Concave	-75mm	1
1"	Bi-Concave	-50mm	1
1"	Bi-Concave	-25mm	1

Hints: Remember to use larger lenses at the input of your system (the side facing the world), and smaller lenses close to the eye. Also, before you start designing, determine which element in the ABCD matrix must equal zero for angular magnification to occur.

**PL1.3 Extra Credit (not required)**: Design a zoom-optic system with the components listed above. (With a zoom-optic system, you can adjust the magnification of your telescope by moving a central lens while keeping the outer most lenses at *approximately* the same location – basically the same setup as is used in an SLR camera zoom lens.)

---

“Two brothers bought a cattle ranch and named it "Focus." When their father asked why they chose that name, they replied: "It's the place where the sons raise meat."

-- Prof. W. B. Pietenpol, Physics Department, University of Colorado, Boulder, Colorado

---

# IN-LAB EXERCISES

## 1.1 Linear Polarizer

### Pre-Lab Question

Explain using electric dipole arguments and with the aid of diagrams why the light reflected from the floor wax is predominantly polarized perpendicular to the plane of incidence.

### In-Lab Assignment

As was discussed in lecture, light reflected off the waxed floor tiles is predominantly polarized perpendicular to the plane of incidence. Holding the given sheet of linear polarizer close to your eye, use such reflections to determine the axis of the polarizer. Convince the T.A. that you have identified the polarizer axis.

T.A. or L.A.  
Signature

## 1.2 Two-Polarizer Diagnostic System

The setup you are given consists of a randomly polarized (unpolarized) laser beam of intensity  $I_0$ , a pair of polarizers, and a screen arranged as illustrated in Figure 1. The second polarizer is often called an analyzer. Its axis is generally oriented perpendicular (crossed) or parallel to the axis of the first polarizer.




Figure 1 – Diagnostic two-polarizer system

### In-Lab Assignment

With the first and second polarizers crossed in the above system, place a third sheet of linear polarizer in the system between the two stationary crossed sheets and observe the intensity of the transmitted light on the screen as you rotate the third polarizer. Make sure to record the orientation of every element in the diagnostic system in your lab notebook.

Now place a photodetector at D and record below the intensity on the screen as a function of the angle of rotation of the third polarizer. About 10 points between (and including) two adjacent intensity maxima should be sufficient. Record this chart into your lab notebook.

Angle	Relative Intensity

**Post-lab Questions**

- (a) Explain why the intensity on the screen changes as the central polarizer is rotated?
- (b) Derive the form of the theoretically expected intensity at the detector, and compare it with your measurement (you may use Jones matrices for this, or use basic trigonometry). Explain the discrepancy, if any.

**1.3 Quarter and Half Wave Plates**

**In-Lab Assignment**

- (a) You are given one quarter-wave plate and one half-wave plate (labeled A and B) but it is not known which plate is which. Use the system in Figure 1 to: (i) find the axes of the quarter-wave plate, (ii) find the axes of the half-wave plate, and (iii) distinguish one plate from the other. Convince the T.A. that you have indeed found the axes and distinguished the two wave plates.

T.A. or L.A.  
Signature

Describe the experiments you performed on the two wave plates. Tabulate the results, and explain how you arrived at your conclusions. Make sure to record the orientation of every element in the diagnostic system in your lab notebooks.

**Post-lab Questions**

A half-wave plate can be used to convert vertically polarized light into horizontally polarized light, and vice versa. Another useful application of a half-wave plate is as a linear polarization rotator. Indeed, if linear polarized light is incident at normal incidence on a half-wave plate and the plate is rotated through an angle  $\theta$  about an axis parallel with the propagation direction, the polarization of the transmitted beam is rotate by  $2\theta$ .

- (a) In your lab notebook, draw diagrams and explain how a half-wave plate converts vertically polarized light to horizontally polarized light,
- (b) In your lab notebook, draw diagrams and explain how and why a half-wave plate achieves polarization rotation.

Note: We are not looking for detailed calculations, but instead a detailed written description using diagrams as visual aids.

## 1.4 Magic Sheet


### In-Lab Assignment

You are given a sheet of "magic" material, which transmits light when folded in one direction, but blocks light when folded in the other direction. Observe this for yourself. Using only linear polarizers, devise nondestructive experiments to determine the possible optical elements that constitute the "magic" material. Record your data in tabular form in your lab notebook, and explain with the help of diagrams why light is blocked or passed in each experiment you perform. Show clearly how you arrive at your conclusion about the magic material.

## 1.5 Dielectric Reflection and the Brewster Angle

### In-Lab Assignment

The Brewster angle is defined as the angle of incidence for which the reflected light intensity from a dielectric boundary goes to zero for TM polarization. This effect can be demonstrated by the experiment shown schematically below.



**Figure 2. - Experimental setup to demonstrate Brewster's angle.**


Use a randomly polarized laser beam for this experiment. Align the glass slide so that its axis of rotation is perpendicular to the laser beam. This can be accomplished by adjusting the laser or the glass slide so that at normal incidence, the reflected beam goes back close to or into the laser cavity.

- (a) Where is the plane of incidence in your experiment?
- (b) With the linear polarizer aligned with its axis horizontal, and the laser beam centered on the axis of rotation of the glass slide, rotate the glass slide and observe the intensity of the reflected beam on the screen. Rotate the polarizer sheet by  $90^\circ$  and repeat the experiment. Copy the table below into your lab notebook. Record your observations of the **reflected beam** for each of the two polarization cases in your lab notebook.



## 1.6 Total Internal Reflection System

You are given a  $45^\circ$ -right-angle glass prism on a rotary table and a HeNe laser set up in the arrangement shown in the Figure 3. Illuminate one of the small faces with the laser beam at normal incidence (reflection off this surface goes back into the laser cavity) and observe the phenomenon of total internal reflection.



**Figure 3 – Setup for observing total internal reflection using  $45^\circ$  - right angle prism.**

### Pre-Lab Questions

- On a sketch of Figure 4, indicate which direction you would have to rotate the prism to get light transmitted through the hypotenuse.
- Derive an expression relating the index of refraction and the rotation angle at which you should first observe light transmitted through the hypotenuse.

### In-Lab Assignment

- Sketch the actual system used for this experiment.

Illuminate one of the small faces with the laser beam at normal incidence (reflection off this surface goes back into the laser cavity) and observe the phenomenon of total internal reflection. Slowly rotate the prism until the transmitted beam just grazes the hypotenuse of the prism.

- Measure the angle of rotation, and use it to calculate the refractive index of the glass block. You should have come up with an equation to do this in the pre-lab exercises. Use your derived equation, while in the Lab to perform a reality check on the index of refraction.

## 1.7 Imaging Properties of a Convex Lens

In this exercise, we will explore the imaging properties of a simple symmetric biconvex lens. The setup we will use is shown in Figure 4. It consists of a white light source a transmission object (slide), a biconvex lens, and a screen, all mounted along an optical rail.




Figure 4. - Setup to explore the imaging properties of a convex lens.

### In-Lab Assignment

Using the setup in Figure 4, and keeping the light source (slide) and object (lamp) fixed, select five different object distances,  $O$ , (distance from the lens to the object) by moving the lens. For each object distance, measure: (1) the corresponding image distance,  $I$ , for sharply defined image on the screen, and (2) the size of the image.

### Post-Lab Questions

- Make a plot of  $1/I$  versus  $1/O$  and from your plot determine the focal length of the lens (use all your data points – if you are not sure how to fit the data, please ask). Include the raw data, the plot and your calculation of the focal length in your report.
- If the refractive index of the lens is 1.52, what is the radius of curvature of its convex surfaces?
- Why is a sharp image not seen when the object distance is less than the focal length of the lens?
- If two lenses of focal length  $F_1$  and  $F_2$  are abutted, what is the focal length of the combination? Show your calculation.
- How would you go about measuring the focal length of a biconcave lens?

## 1.8 Telescope Fabrication (optional - for extra credit)

Build the telescope you designed in the pre-lab exercise. Comment on the performance of your telescope. Did your telescope operate as you expected? Why or why not? Explain the operation of your telescope to the TA.

T.A. or L.A.  
Signature

## APPENDIX I

As mentioned earlier part of the Lab 1 write-up will be an oral presentation. You have two choices for topics. For the topic you choose you should discuss the relevant experiments as performed in lab and answer the related post lab questions.

### (1) Wave plates and polarization (Experiments 1.1 to 1.3)

This topic deals with polarization and how waveplates can be used to change polarization. Your presentation should discuss what polarization is and how waveplates change the polarization of light passing through them. You should present the experiments, the motivation behind the experiments (what each experiment was supposed to show), and how your group determined the composition of the "Magic Sheet".

### (2) The Fresnel Reflection Coefficients (Experiment 1.5)

This topic deals with the Fresnel reflection coefficients which determine how the intensity of a reflected beam changes with incidence angle. Your presentation should discuss the Fresnel reflection coefficients in general and how they can be used to determine the index of a material given a polarized source. You should be able to discuss why the coefficients are different for TE and TM polarizations, and explain what is happening at Brewsters angle.

You should prepare a ten minute presentation with no more than 8 slides and be ready for 5 min of discussion on the topic you choose to present. Keep in mind that 8 is an upper limit, and you should really have 4-6 content slides. Your presentation should cover the post-lab questions from the lab assignment for the topic you choose, and should be aimed as if you were presenting to the other 6.161 students (who don't necessarily understand the lab material).

This cannot be said enough, but practice your presentation. It does make a difference.

The post lab questions for the rest of the experiments should be written up in your lab notebook, and I'll check them after lab next week.

## APENDIX II

Student:  
 Date:            Lab:  
 Grade:  
 CI Instructor: Don Unger

### 6.161 CI-M, Oral Presentation Rubric, Fall 2011

	Okay	Good	Excellent	Comments
<b>Content</b>				
Introduces self; purpose of presentation is clearly stated.	8	10	12	Category Max 36 points
Presentation has a beginning/overview, middle, and end.				
Points are supported with clear examples.				
<b>Style and Delivery</b>				
Maintains eye contact with audience.	2	3	4	Category Max 28 points
Avoids filler words (e.g., <i>um</i> , <i>like</i> , ah)				
Language appropriate, professional				
Modulates voice appropriately (pitch, volume, inflection, tone, pace)				
Body language indicates interest, professionalism, and poise (gestures, fidgets, posture, clothing choice).				
Pointer <sup>1</sup> used effectively				
Honors time limitations.				
<b>Slides</b>				
Readable from back of room	7	8	9	Category Max 36 points
Titles orient and focus audience				
Data presentation: legible, uncrowded, properly labeled, symbols distinct, error bars and plot points present; use of color allows for lighting conditions, projector distortion and weaker eyes				
Slides contain only minimal text				

1. Staff will provide a combined laser pointer, slide advance remote.

V 1.0 3Jul2011 DNSU

-----

Congratulations, you have finished Laboratory #1. Please write down in your notebook how long it took (start to finish) to complete: (a) the pre-lab work, (b) the in-lab portion, and (c) the post-lab writeup.