

DCIF NMR TRAINING SUMMARY

600 MHz

BRUKER DRX "B600"

last edit 20071101

TOPSPIN COMMANDS

Instrument Setup

- | | |
|---|-----------------------|
| 1. Create dataset | edc or new |
| 2. Turn on eject gas | ej |
| 3. Insert sample and turn off eject gas | ij |
| 4. Turn spin gas on at 20 Hz | ro on |
| 5. Set temperature | edte |
| 6. Open lock display | lockdisp |
| 7. Load bestshim file | rsh best |
| 8. Turn lock on | lock & choose solvent |
| 9. Open BSMS display | bsmsdisp |
| 10. Shim sample in BSMS then close BSMS | |
| 11. Gradient Shimming | gradshim |

Experiment Setup and Acquisition

- | | |
|---|----------------|
| 12. Load parameter set (<i>if not loaded in edc window</i>) | rpar |
| 13. Load prosol parameters (p1, pl1, etc) | getprosol |
| 14. Adjust parameters (ns, sw, d1, etc) | ased |
| 15. Set solvent | solvent |
| 16. Setup console | ii |
| 17. Check tuning | wobb |
| 18. Set receiver gain | rga or rgacryo |
| 19. Start experiment | zg |
| 20. Transfer acquired FIDs | tr |

Data Processing

- | | |
|---|------|
| 21. Fourier Transform | ft |
| 22. Autophase spectra | apk |
| 23. Stop experiment and save data | halt |
| 24. Invoke exponential function, FT and phasing | efp |
| 25. Baseline correction | abs |
| 26. Start calibration menu | cal |
| 27. Start peak picking menu | pp |
| 28. Start integration menu | int |
| 29. Print | |
| a. print active screen | prnt |
| b. open plot editor | plot |

Misc

- | | |
|---|---------------------|
| • Increment the experiment number | iexpno |
| • Read dataset | re <expno> <procno> |
| • Run sequential experiments | multizg |
| • Resume acquisition and add to dataset | go |

When Finished

- | | |
|------------------------|--|
| 30. Turn lock off | lock off |
| 31. Turn spin off | ro off |
| 32. Eject sample | ej |
| 33. Turn off eject gas | ij |
| 34. Exit program | exit |
| 35. Log off computer | Right click desktop choose [Logout] |
| 36. Sign log book | |