

DCIF TRAINING GUIDE

KINETICS

BRUKER 400/401/600

(20080919)

As with all of the other DCIF tutorials, the conventions used in this guide are as follows:

- **Boldface** type indicates commands that are typed into the Topspin command line or in a terminal window
- Italic [*Boldface*] type with square brackets indicates a button in the Topspin menu that is to be pushed.
- <**Boldface**> type surrounded by a bracket indicate keyboard strokes.
- **LMB** indicates the Left Mouse Button
- **MMB** indicates the Middle Mouse Button
- **RMB** indicates the Right Mouse Button
- Unless otherwise noted, all commands typed into the input window are followed by an <↵ Enter> keystroke.

1. Insert, lock and shim your sample.
2. Setup the parameter set; tune the probe if necessary.
 - a. In order to measure accurate integrals you should determine the T1 relaxation time for your molecule. Your delay should be 5x longer than T1 for quantitative measurements.
3. Acquire a trial run and note how long the acquisition takes.
4. Type **dcif-kinetics**
 - a. “This program will automatically increment the experiment number for a given number of acquisitions.” Continue (y/n) **y**
 - b. “How many data sets will you acquire (total)?” **#**
 - c. “How long of a delay (in seconds) between experiments?” **#** (seconds)
 - i. This is the delay between when each experiment starts. Make sure this delay is longer than the total acquisition time needed.
 - d. The acquisition should start automatically.
5. When the entire array is complete Topspin should display “ Kinetics Experiment Complete”
6. The resulting experiment numbers will be contiguous. If you setup and start the array of 10 spectra in expno 1, expno 2 through 11 will have your data.
 - a. The command **dpa** will display the time the acquisition started.
 - b. To review each experiment type **re <#>** or type **search**.
7. Use the TopSpin Plot Editor to print stacked spectra.