



Extending Microarray Technology to Study Protein Function

Gavin MacBeath
Bauer Center for Genomics Research
Harvard University

MICROARRAYS: THE POWER OF MULTIPLEXING





MICROARRAYS: THE POWER OF MULTIPLEXING





KEEP THE PROTEINS HAPPY

- hydrophilic surface
- amine-reactive chemistry
- constant hydration


BSA-NHS SLIDES




ALDEHYDE SLIDES




PROTEIN MICROARRAYS



PROTEIN-PROTEIN INTERACTIONS




PROTEIN-PROTEIN INTERACTIONS




0.50 $\mu\text{g/ml}$ BODIPY-FL-IgG
 0.05 $\mu\text{g/ml}$ Cy3-IkB α
 0.50 $\mu\text{g/ml}$ Cy5-FKBP12

PROTEIN-PROTEIN INTERACTIONS




0.50 $\mu\text{g/ml}$ BODIPY-FL-IgG
 0.50 $\mu\text{g/ml}$ Cy5-FKBP12 + 100 nM rapamycin


TITRATION ON GLASS SLIDE




PROTEIN MICROARRAYS




ENZYME-SUBSTRATE INTERACTIONS




IDENTIFYING SUBSTRATES OF PROTEIN KINASES




PROTEIN MICROARRAYS




PROTEIN-SMALL MOLECULE INTERACTIONS



IDENTIFYING THE TARGETS OF SMALL MOLECULES



PROTEIN MICROARRAYS



FUNCTIONAL GENOMICS

subsets of related proteins

- informatics (sequence motifs)
- literature searches

protein domains

- coiled coils
John Newman
- protocadherins
Viara Grantcharova

full-length proteins

- cancer proteins
Joshua LaBaer, Pascal Braun
- c. elegans* development
Viara Grantcharova

COILED COILS



myc/max
heterodimer

J Mol Biol (1998) **281**, 165.

Mediate homo- and hetero-dimerization

Found in:

- structural proteins
- motor proteins
- transcription factors
- membrane fusion proteins

In yeast, MULTICOIL predicts:


- ~300 proteins with 2-stranded coiled-coils
- ~250 proteins with 3-stranded coiled coils

tropomyosin
homodimer

J Mol Biol (1986) **192**, 111.




COILED-COIL SCREEN: YEAST TWO-HYBRID




Newman, Wolf & Kim (2000) *PNAS* **97**, 13203

COILED-COIL INTERACTIONS: PROTEIN MICROARRAYS




John Newman

PILOT COILED-COIL SCREEN: PROTEIN MICROARRAYS




John Newman


PROTEIN MICROARRAYS




SCREENING METHOD 2




SCREENING METHOD 2




SCREENING METHOD 2




SCREENING METHOD 2



SCREENING METHOD 2



SCREENING METHOD 2



THE MICROARRAY WORLD



Glass slides

- 2.5 cm x 7.5 cm
- >10,000 spots per plate
- spacing varies


THE HTS WORLD




Microtiter plates

- 8.5 cm x 12.5 cm
- 96, 384, or 1536 wells per plate
- 9 mm, 4.5 mm, or 2.25 mm spacing

MICROARRAYS IN WELLS OF PLATES




MULTIPLEXED SCREENING




- bottomless 384-well plate
- 64 wells per slide x 4 slides per plate = 256 wells per plate
- 256 wells per plate x 100 proteins per well = 25,600 assays per plate

PROTEIN MICROARRAYS




ANTIBODY ARRAYS: SANDWICH ELISA




ANTIBODY ARRAYS: MICRO-SANDWICH ASSAY

Flourescent 2' Antibodies




MICRO-SANDWICH DETECTION TUMOR CELL PROFILING




ANTIBODY ARRAYS: MICRO-SANDWICH ASSAY

Ratiometric analyses of protein modification




PHOSPHORYLATION AND REGULATION OF GROWTH FACTOR RECEPTORS



- EGFR is phosphorylated in response to EGF
- ErbB2 is phosphorylated by EGFR
- EGFR is down-regulated
- TfR remains unchanged

RATIOMETRIC ANALYSES OF EXPRESSION AND PHOSPHORYLATION





Capture and detection antibodies against different epitopes

Blue - Alexa488; Green - Cy3; Red - Cy5

RATIOMETRIC PROFILING IN TUMOR CELL LINES



EFFECT OF EGF



- RED = growth factor receptor expression
- GREEN = growth factor receptor phosphorylation
- BLUE = transferrin receptor expression
- YELLOW = ratio expression/phosphorylation


RATIOMETRIC ANALYSES OF EXPRESSION AND PHOSPHORYLATION

EGF time course on SK-BR-3 cells




RATIOMETRIC ANALYSES OF EXPRESSION AND PHOSPHORYLATION

EGF time course on SK-BR-3 cells



PROTEIN MICROARRAYS



SUMMARY

- rapid discovery of protein interactions
- screening of small molecules *in solution*
- screening for *specificity*
- system-wide analysis at the *protein level*

ACKNOWLEDGEMENTS

Protein Microarrays

Prof. Peter Kim (Whitehead Institute)
Dr. John Newman (Whitehead Institute)

Dr. Viara Grantcharova
Lioudmila Zaslavskaja

Antibody Microarrays

Prof. Peter Sorger (MIT)
Dr. Ulrik Nielsen (MIT)
Dr. Michael Cardone (MIT)

Dr. Raghida Bu-Khalid

Harvard Center for Genomics Research
DARPA