

Instructions for SecureFX 6.0.x

<https://ist.mit.edu/services/software>

How to backup your files

last revised 1/28/2010

Setup:

Select **File** → **Connect** to open the Connect pop-up window.

Click the New Session icon (third one in from the left) to open the New Session Wizard.

Enter the following parameters:

SecureCRT protocol:	SSH2
Host Name	IP address (18.60.XX.XX)
Port:	22
Firewall:	None
User name	<i>Leave blank</i>
SecureFX protocol:	SFTP
Session name:	<i>Whatever you wish to name the session</i>
Description:	<i>Again, whatever you wish</i>

Click Finish

Connecting to a Spectrometer Computer

Select **File** → **Connect** to open Connect pop-up window.

Highlight the desired computer (session) and click Connect

Enter your User ID and password for that computer.

Windows for the remote and local hosts will display the associated file tree.

File Transfer

Method A: Manual Download

Highlight the folder or file to download.

Select **File** → **Manual Download**

Enter the local host location for the file.

Enter Binary as the file transfer type.

Click OK

Method B: Right Click Drag and Drop

Highlight the remote folder or file to download.

Hold the LEFT mouse button down and drag the file to local host window.

Release the LEFT mouse button. If prompted, select Download Here as Binary.